



Concentrate for disinfectant cleaning of medical devices and other surfaces

terralin[®] protect

Our Plus

- aldehyde free
- excellent cleaning performance
- discreet, pleasant and fresh odour

Application areas

Disinfection and cleaning of medical devices and other surfaces:

- In all hospital areas (e.g. wards, ICU, theatres)
- Particularly for odour-sensitive areas (e.g. neo-natal intensive care units)
- On sensitive materials (e.g. acrylic glass)

Instructions for use

Prepare use-solution in concentration required. Make sure to wet surfaces completely and keep them wet for the whole exposure time. Exact dosages can be achieved with the use of schülke dosage aids. Wipe surfaces with a damp cloth. Always aim for complete surface coverage. Rinse equipment carefully with water. In order to avoid loss of microbiological effectiveness do not mix this product with any other cleaners. Therefore, please seek advice from schülke before attempting such mixtures. Not for final disinfection of semicritical and critical medical devices!

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>EN1276, EN13697, in accordance with VAH - dirty conditions</i>	0.5 % (5 ml/l)	60 min.
bactericidal <i>in accordance with VAH</i>	2 % (20 ml/l)	15 min.
tuberculocidal <i>in a suspension test</i>	0.5 % (5 ml/l)	2 h

Efficacy	Concentration	Contact time
yeastcidal <i>EN1650, EN13697, in accordance with VAH</i>	0.5 % (5 ml/l)	60 min.
virucidal against enveloped viruses <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	0.5 % (5 ml/l)	60 min.
Adenovirus <i>in accordance with DVV</i>	2 % (20 ml/l)	4 h
Norovirus <i>- clean conditions</i>	1 % (10 ml/l)	2 h
Norovirus <i>- dirty conditions</i>	2 % (20 ml/l)	2 h
Polyoma SV40	2 % (20 ml/l)	2 h
Rotavirus <i>in accordance with DVV</i>	0.25 % (2.5 ml/l)	5 min.

Application area	Concentration	Contact time
bathtub disinfection	ready-to-use	5 min.

Certificates

- IHO listed
- VAH certificate



Product data

Composition:

100 g solution contains the following active substances:
22 g Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-C16)), 17 g 2-Phenoxyethanol, 0.9 g Amines, N-C12-C14 (even numbered)-alkyltrimethylenedi-, reaction products with chloroacetic acid)

Labelling according to Regulation (EC) No. 648/2004:

5 - 15 % non-ionic surfactants, perfumes

Chemical-physical data

Color	green
Density	ca. 1,01 g/cm ³ / 20 °C
Flash point	48 °C / Method : DIN 51755 Part 1
Form	liquid
pH	ca. 8,6 / 20 °C
Viscosity, dynamic	ca. 21 mPa*s / 20 °C / Method : ISO 3219

Special advice

Use disinfectants safely. Always read the label and product information before use.

Information for order

Item	Delivery form	Item no.
terralin [®] protect 2 l bottle	5/Carton	on request
terralin [®] protect 200 l	1/drum(s)	on request
terralin [®] protect 500 l	1/Container	on request
terralin [®] protect 20 ml bag (500x20ml)	1/Carton	on request
terralin [®] protect 5 l	1/Canister	on request
terralin [®] protect 20 l	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Can key for 20 + 25 l	135901
Can key for 5 + 10 l	135810
Decanter funnel	117901
measuring cup 500 ml	136101
Opener all-in-one (1l eurobottle/can/container)	70000804
schülke dosing feeder 20 l / 25 l (20ml / stroke)	182501
schülke dosing feeder 5 l / 10 l (20ml / stroke)	117101
schülke tap for 5 l / 10 l can	135501
tap for 20 / 25 lt can	126101

Related Products

- s&m care additive
- s&m cleaning additive
- schülke wipes
- schülke wipes safe&easy

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com