

Document Number: V200QARA-SWI-01-A

Revision Level: 02

TITLE: Technical Data Sheet

Page 1 of 13

# BD PosiFlush<sup>™</sup> SP & XS Syringe pre-filled with a sterile 0.9% Sodium Chloride Solution, Sterile

BD Switzerland Sàrl Terre Bonne Park – A4 Route de Crassier 17 1262 Eysins, Switzerland bd.com

TDS number: V201-015 - Rev. 02

2021-February

## 1. General Information

# 1.1 <u>Intended use</u>

The BD PosiFlush™ SP and XS syringe are intended for immediate use in maintaining patency of vascular access devices.

BD PosiFlush™ SP and XS syringe are not intended for dry product reconstitution, for medication dilution, or where intravenous therapy with sodium chloride is indicated.

# 1.2 **General description**

The syringe consists of a ready to use syringe, which is pre-filled with an injectable sterile, non-pyrogenic and isotonic 0.9% sodium chloride solution. The syringe is not to be resterilized or re-used.

The BD PosiFlush<sup>TM</sup> XS (externally Sterile) is presented with individual protection for maintaining sterility and can be used on a sterile field. BD PosiFlush<sup>TM</sup> SP cannot be used on a sterile field.

All sizes of the pre-filled syringes (3, 5 and 10 ml) share a common barrel diameter. The tip cap is the same for each type of syringe. The plunger stopper is identically assembled on each size of syringes (3, 5, 10 ml).





Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 2 of 13

	BD Catalog Number	BD Product Description	Color Code
	306573	3ml BD PosiFlush™ SP	Orange
BD	306574	5ml BD PosiFlush™ SP	Orange
PosiFlush™	306575	10 ml BD PosiFlush™ SP	Orange
SP, sterile	306583	3ml BD PosiFlush™ EMA SP	Orange
Fluid Path	306584	84 5ml BD PosiFlush™ EMA SP	
	306585	10 ml BD PosiFlush™ EMA SP	Orange
DD.	306570	3ml BD PosiFlush™ XS	Blue
BD PosiFlush™	306571	5ml BD PosiFlush™ XS	Blue
XS,	306572	10 ml BD PosiFlush™ XS	Blue
externally	306580	3ml BD PosiFlush™ EMA XS	Blue
sterile	306581	5ml BD PosiFlush™ EMA XS	Blue
Sterne	306582	10 ml BD PosiFlush™ EMA XS	Blue

Note: Please check BD catalog number availability in your country.

The BD Product Description can slightly differ from the Declaration of Conformity; please always refer to the BD Catalog Number.

# **Further features:**

N/A

# 1.3 <u>Certification</u>

BD	BD Legal Manufacturer	CE Certificate Number And Notified Body Brief Name  CE Certificate BD Manufacturing Site (Country of Origin) and ISO 13485 Certification		EC
Catalog	and			Representative
Number	ISO 13485 Certification			(if applicable)
306573 306574 306583 306584	Address: Becton Dickinson and Company 1 Becton Drive Franklin Lakes New Jersey 07417 United States  ISO 13485 Certificate No.: MD19.2305	CE certified with NSAI (0050) Certificate No.: 252.780 AND CE certified with GMED (0459) Certificate No.: 16730	Address: BD Medical Surgical 2153 12th Avenue Columbus NE 68601 United States  ISO 13485 Certificate No.: MD19.2143  OR  Becton Dickinson and Company 12 Av Mequinenza 22520 Fraga, Hesca, Spain ISO 13485 Certificate No.: 2015 05 0047 EN	Becton Dickinson and Company Donore Road Drogheda Co Louth, Ireland





Document Number: V200QARA-SWI-01-A

TITLE: Technical Data Sheet

Revision Level: 02 Page 3 of 13

BD Catalog Number	BD Legal Manufacturer and ISO 13485 Certification	and And (Country of Origin)		EC Representative (if applicable)
306575 306585	Address: Becton Dickinson and Company 1 Becton Drive Franklin Lakes New Jersey 07417 United States  ISO 13485 Certificate No.: MD19.2305	CE certified with NSAI (0050) Certificate No.: 252.780 AND CE certified with GMED (0459) Certificate No.: 16730	Address: BD Medical Surgical 2153 12th Avenue Columbus NE 68601 United States  ISO 13485 Certificate No.: MD19.2143  OR  Becton Dickinson and Company 12 Av Mequinenza 22520 Fraga, Hesca, Spain ISO 13485 Certificate No.: 2015 05 0047 EN  OR  Address: Becton Dickinson and Company Donore Road Drogheda Co Louth, Ireland  ISO 13485 Certificate No.: MD19.1609	Becton Dickinson and Company Donore Road Drogheda Co Louth, Ireland
306570 306571 306572 306580 306581 306582	Address: Becton Dickinson and Company 1 Becton Drive Franklin Lakes New Jersey 07417 United States  ISO 13485 Certificate No.: MD19.2305	CE certified with NSAI (0050) Certificate No.: 252.780 <b>AND</b> CE certified with GMED (0459) Certificate No.: 14879	Address: Becton Dickinson and Company Donore Road Drogheda Co Louth, Ireland  ISO 13485 Certificate No.: MD19.1609  Becton Dickinson and Company Donore Road Drogheda Co Louth, Ireland	



Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 4 of 13

## 1.4 Materials

Component	Material
Saline Solution	Sodium Chloride 0.9%
Barrel	Polypropylene
Tip Cap	Polypropylene + Colorant
Stopper	Elastomer
Plunger Rod	Polypropylene
Lubricant	Silicon oil

# 1.5 Materials of concern

Materials of concern are chemicals or substances that have been identified as having the potential to cause long term effects on humans or the environment.

Material	Comment		
Phthalates	Based on our ongoing data collection efforts and/or information received from our suppliers as per 10 June 2020, BD has not identified any  1,2-Benzendicarboxylic acid, dihexyl ester (branched & linear) (CAS#68515-50-4),  1,2-Benzendicarboxylic acid, di-C6-8-branched alkyl esters (CAS#71888-89-6),  1,2-Benzendicarboxylic acid, di-C7-11-branched and linear alkyl esters (CAS#68515-42-4),  1,2-Benzendicarboxylic acid, di-C6-10 alkyl esters (CAS#68515-51-5),  1,2-Benzendicarboxylic acid, mixed decyl, hexyl, and octyl diesters (CAS#68648-93-1),  Benzyl butyl phthalate (BBP) (CAS# 85-68-7),  Bis (2-ethylhexyl) phthalate (DEHP) (CAS# 117-81-7),  Bis (2-methoxyethyl) phthalate (DMEP) (CAS# 117-82-8),  Di-n-hexyl phthalate (DnHP) (CAS# 84-75-3),  Dibutyl phthalate (DBP) (CAS# 84-69-5),  Diisopentylphthalate (DIBP) (CAS# 84-69-5),  Diisopentylphthalate (DIPP) (CAS# 84-69-5),  Dipentyl phthalate (DPP) (CAS# 84-69-9), or  Dicyclohexyl phthalate (CAS# 776297-69-9), or  Dicyclohexyl phthalate (DCHP) (CAS# 84-61-7)  in the articles and packaging with the Product Numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w).		
Latex	Based on our ongoing data collection efforts and/or information received from our suppliers as per 10 June 2020, the articles with the Product Numbers above are not formulated with natural rubber latex.		
Bisphenol A	Based on our ongoing data collection efforts and/or information received from our suppliers as per 10 June 2020, BD has not identified any  • 4,4'-isopropylidenediphenol (BPA) (CAS# 80-05-7) in the articles and packaging with the Product Numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w). It is not a building block of any of the raw materials utilized and is not intentionally added. BD has not done any testing to evaluate levels of this chemical in these products.		
Substances of animal origin BSE/TSE	The raw materials used in the manufacture of these devices do not contain any animal tissue but may contain very small amounts of chemicals that are derived from animal-derived raw materials. These products are manufactured using polymer resins which may contain very small amounts of stearic acids and related substances derived from tallow derivatives. Our resin suppliers have confirmed that these chemicals have been produced		



Document Number: V200QARA-SWI-01-A

Revision Level: 02

TITLE: Technical Data Sheet

Page 5 of 13

	with multiple cycles of conditions at least as rigorous as those specified in Annex C.5 of EN ISO 22442-1:2015 and Section 6 of EMA 410/01 Rev. 3. Therefore, these raw materials meet or exceed the requirements of EN ISO 22442-1 and EMA 410/01 Rev. 3. Based on this information, these products are considered not to present any risk with respect to TSE/BSE or other animal-borne diseases. Furthermore, as recognized by MEDDEV 2.4/1, such derivative chemicals produced in accordance with the aforementioned standards and guidelines are considered irrelevant when determining the classification of a medical device (per MDD 93/42/EEC and EU No 722/2012).
Polyvinyl chloride (PVC)	The medical devices and packaging referenced above have not been designed nor intentionally manufactured with any additives or raw materials that contain PVC. No PVC is intentionally added to these medical devices.

## 1.6 REACH information

Based on our ongoing data collection efforts and/or information received from our suppliers as per 10 June 2020, BD has not identified any chemicals in the articles and packaging with the Product Numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w), which have been listed as SVHC and included in the "Candidate List" published by the European Chemical Agency (ECHA) on 16 January 2020 according to Art. 59 (1,10) of the Regulation (EC) N° 1907/2006 (REACH).

# 1.7 **Biocompatibility**

BD Medical products comply with the requirements of the standard for toxicity, pyrogenicity and biocompatibility of medical devices, ISO 10993 series - Biological Evaluation of Medical Devices.

# 1.8 Sterilization method

The steam sterilization method guarantees a Sterility Assurance Level of 10-6 for the finished product, including fluid path, solution and outer surface of the syringe.

For the BD PosiFlush™ SP and XS product, the external syringe is also sterile to a Sterility Assurance Level of 10-6. The validation and controls are performed according to EN ISO 17665-1 requirements and EN 556-1.

#### 1.9 Shelf life and storage conditions

The BD PosiFlush™ SP and XS shelf life have been assessed by stability studies in order to verify the functionality, physico-chemical and microbial properties over time. BD PosiFlush™ SP and XS have a shelf life of 3 years.

The resulting recommendations in terms of transportation and storage of the sodium chloride pre-filled syringes are:

- They need to be shipped in order to maintain the conditions of transportation of 0-40°C (32-104°F) targeting 25°C
- They need to be stored under specific conditions: the temperatures need to be maintained at 15°C 25°C (59° 77°F) with excursions permitted to 30°C (86°F)
- They must not freeze





Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 6 of 13

#### 1.10 Standards

As per extract from the Declaration of Conformity DoC PosiFlush SP and XS rev B linked to CE certificate number 252.780:

CL Certificate Humber	232.700.
Harmonized Standards	
EN 556-1:2001	Sterilization of Medical Devices - Requirements for Medical Devices to Be Designated "Sterile" - Part 1: Requirements for Terminally Sterilized Medical Devices
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN 62366:2008	Medical Devices – Application of usability engineering to medical devices
EN 22442:2007	Medical devices utilizing animal tissues and their derivatives
EN ISO 15223-1:2016	Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General Requirements
EN ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-4:2009	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-11:2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
EN ISO 10993-17:2009	Biological evaluation of medical devices - Part 17: Establishment of allowable limits for leachable substances
EN ISO 10993-18:2009	Biological evaluation of medical devices - Part 18: Chemical characterization of materials
EN ISO 13485:2016	Medical devices – Quality management systems – Requirements for regulatory purposes
EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
EN ISO 17665-1:2006	Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
EN ISO 1707:1996	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Lock fittings
EN 20594-1:1993	Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment – Part 1: General requirements
EN ISO 11737-1:2006	Sterilization of medical devices – Microbial methods- Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2009	Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
Non-Harmonized Stand	dards
BS EN 540:1993	Clinical investigation of medical devices for human subjects
ICH Q1A Guidance	Guideline for stability testing of next drug substances and products
ISO 14644:2000	Clean rooms and associated controlled environments
USP<42>	Version 42
ISO 10993-10:2010	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization

## Note:

The above standards reflect the status at the time of drafting this document. More information or updates are available on request in the Declaration of Conformity.



Document Number: V200QARA-SWI-01-A

Revision Level: 02

TITLE: Technical Data Sheet

Page 7 of 13

#### 1.11 Classification

The device is a pre-filled syringe with a solution, but as it is to be used for a mechanical action, it is a Class III medical device under rule 13 of Annex IX of the Directive 93/42/EEC.

#### 1.12 GMDN code

According to ISO 15225 (Medical devices - Quality management - Medical device nomenclature data structure), BD PosiFlush™ SP and XS are referenced as follows:

GMDN Code: 64786

GMDN Term: Vascular Catheter/Cannula flush solution, non-anticoagulant, non-antimicrobial

# 1.13 Manufacturing practices

The entire manufacturing and testing processes are following the Manufacturing Practices as specified below:

- Incoming raw materials are verified via material inspection and testing and our suppliers are approved via our vendor management system.
- In addition to the automatic on-line inspections, in-process inspections are performed in addition to final product testing to ensure compliance with approved specifications.
- The manufacturing and testing details of each batch of product are recorded on a batch record which is retained in accordance with our document control procedures.
- BD operates a system of internal and external audits to maintain compliance.
- BD confirms that it will continue to adhere to relevant international standards in designing and manufacturing its products.
- BD reserves the right to use the internal change control procedure to change raw material suppliers and production process.

#### 1.14 Other information

- (Material) Safety Data Sheets are not required for this product.
- Certificate of Food Contact (Commission Regulation EU 1183/2012 on "plastic materials and articles intended for contact with food" and Directive 2002/72/CE (as amended) "relating to plastic materials and articles intended to come into contact with foodstuffs") is not required as BD products are used for general purpose injection and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.
- Good Manufacturing Practices as defined by the FDA Pharmaceutical is not applicable for Medical Devices.
- The EC representative is BD Drogheda, Ireland



Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 8 of 13

2. Packaging

# 2.1 Packaging configuration

	BD Catalog Number	BD Product Description	Primary Packaging (Qty)	Shelf Box (Qty)	Shipping Case (Qty)	IFU Insert N/A / Yes / No*
BD	306573	3ml BD PosiFlush™ SP	1	30	480	Yes
PosiFlush	306574	5ml BD PosiFlush™ SP	1	30	480	Yes
™ SP,	306575	10 ml BD PosiFlush™ SP	1	30	480	Yes
sterile	306583	3ml BD PosiFlush™ EMA SP	1	30	480	Yes
Fluid Path	306584	5ml BD PosiFlush™ EMA SP	1	30	480	Yes
Fiulu Patii	306585	10 ml BD PosiFlush™ EMA SP	1	30	480	Yes
	306570	3ml BD PosiFlush™ XS	1	30	240	Yes
BD Boo!Elwah	306571	5ml BD PosiFlush™ XS	1	30	240	Yes
PosiFlush	306572	10 ml BD PosiFlush™ XS	1	30	240	Yes
™ XS,	306580	3ml BD PosiFlush™ EMA XS	1	30	240	Yes
externally sterile	306581	5ml BD PosiFlush™ EMA XS	1	30	240	Yes
Sterne	306582	10 ml BD PosiFlush™ EMA XS	1	30	240	Yes

<sup>\*&</sup>quot;No": IFU may be available but not as an insert.

# 2.2 <u>Packaging material</u>

Component	Material
Unit Pack (Flow wrap) PosiFlush SP	Polypropylene
Unit Pack (Blister Pack) PosiFlush XS	Top Web: Steam paper Bottom Web: Polypropylene
Shelf Box	Chipboard carton
Shipping Case	Corrugated carton



Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 9 of 13

### 2.3 <u>Examples of labeling</u>

Labels: According to European Medical Device directive, labels are multilingual.

Primary Packaging Label extracted from document DG2042 related to reference 306573 (made in Fraga):



Unit pack (flow wrap) extracted from document DGW1511 related to reference 306573 (made in Fraga):

Do not use on a sterile field.
Ne pas utiliser sur champ stérile.
Nicht in einem sterilen Feld verwenden.
No utilizar sobre un campo estéril.
Não utilizar num campo estéril.
Non utilizzare su un campo sterile.
Να μην χρησιμοποιείται μέσα σε στείρο πεδίο.
Niet gebruiken op een steriel veld.
Får ej användas på ett sterilt område.
Ei saa laittaa steriilille alueelle.
Må ikke brukes på et sterilt område.
Må ikke anvendes på et sterilt felt.



Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 10 of 13

Shelf Box extracted from document 10000303830 related to reference 306573 (made in Fraga):



Shelf Box label extracted from document 10000334556 related to reference 306573 (made in Fraga):






Document Number: V200QARA-SWI-01-A

TITLE: Technical Data Sheet

Revision Level: 02 Page 11 of 13

Shipping Case extracted from document 10000303831 related to reference 306573 (made in Fraga):




Case Label extracted from document 10000334557 related to reference 306573 (made in Fraga):





Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 12 of 13

IFU insert (English part only) extracted from document 10000223600 related to reference 306573 (made in Fraga):







Document Number: V200QARA-SWI-01-A Revision Level: 02
TITLE: Technical Data Sheet Page 13 of 13

ex. 7. EXPLICACIÓN DE LOS SIMBOLOS DE 7. PERLOXÃO DOS SIMBOLOS EL 7. PERCAZIONE DEI SIMBOLI
el 7. TELEPHERI TICH SYMBOLAR EL 7. VERRALINO AN DEI SYMBOLER 82. 7. FORTALINO AN SYMBOLER 82. 7. FORTALINO AND SYMBOLER 82. 7. FOR

REVISION	CHANGE SUMMARY
01	Initial release according to new template
02	Update of 1.5: Material of concern Added details in 1.9: Shelf life and storage conditions Update of 1.10: Standards Update of 1.12: GMDN code Update of 2.3: Examples of labelling