

# MICRO-TOUCH<sup>®</sup> Nitrile EP<sup>™</sup>

P e r f o r m a n c e   S e r i e s

## Nitrile, powder-free examination glove with an extended cuff

This non-sterile glove has an extended cuff provides extra protection when healthcare providers are exposed to chemotherapy agents, bodily fluids or blood-borne pathogens in high-risk clinical settings. It provides excellent chemical and puncture resistance and is textured for better grip in wet and dry procedures.

### KEY FEATURES AND BENEFITS:

- Extended 12" cuff provides added protection over the wrist and forearm
- Effective barrier protection against punctures and a wide range of chemicals
- Protects from latex Type I allergy
- 200-count packaging reduces waste
- Tested for use with chemotherapy drugs\*

\* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

\*\* This product is only available in Europe, the Middle East and Africa (EMEA).

### Recommended For

- Tested for use with chemotherapy drugs
- Extra protection over the wrist and arm
- Protection from Type I latex allergy in HCW's or patients
- Sustained, high-risk examination procedures
- Environmentally friendly facilities

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> Nitrile EP<sup>™</sup>

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	White	Cuff Length	Extended
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.07 / 2.8
Finger Thickness (mm/mil)	0.10 / 3.9	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	265 / 10.4	Product Segmentation	Performance Series

## ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	700324	700323	700322	700325	700321

## PACKAGING AND STORAGE

Packaging	200 gloves per dispenser 10 dispensers per case 2000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, ASTM D6978, Category III, EN 374:2003, EN 374:2003, EN 420:2003 + A1:2009, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 421:2010, ISO 10993-10:2002, ISO 13485, ISO 9001

## Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com