



## MICRO-TOUCH® NITRILE

### A NEW WHITE NITRILE GLOVE: OPTIMISED ALLERGY PREVENTION AND COMFORT WITH ECONOMICAL PACKAGING

Today, healthcare workers are highly aware of the risk of Type I allergies resulting from exposure to latex. Made of 100 % nitrile, Ansell's MICRO-TOUCH® Nitrile gloves contain no natural rubber latex and are powder-free.

They provide an effective barrier protection against punctures and a wide range of chemicals, while the white colour reveals soiling and contamination.

MICRO-TOUCH® Nitrile provides latex-like comfort. 'Thermo-elastic' material that adapts to the shape of the hand ensures even greater comfort, while an inner polymer coating makes donning easy and the absence of powder considerably reduces the risk of skin abrasion.

The micro-textured finger surfaces ensure a good grip on instruments. The combination of easy donning, supreme comfort and sure grip make MICRO-TOUCH® Nitrile the optimal choice of many healthcare professionals.

Moreover, the 150-count packaging reduces waste and saves on costs.

This examination glove, CE Marked as Class I Medical Device, is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.

### PRODUCT DESCRIPTION

Material	Nitrile
Colour	White
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Polycrylate coated

### SIZE / PRODUCT CODE

<b>x-small</b>	700151	<b>small</b>	700152
<b>medium</b>	700153	<b>large</b>	700154
<b>x-large</b>	700155		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

### Nitrile Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III  
 Registering Authority: Centexbel (0493): PPE  
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at [www.ansell.eu](http://www.ansell.eu).

### Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J  
 Boulevard International 55, 1070 Brussels, Belgium  
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03  
<http://www.ansell.eu> • E-mail [info@ansell.eu](mailto:info@ansell.eu)

