

# MICRO-TOUCH<sup>®</sup> Nitrile

Essential Series

## Nitrile, powder free examination glove with exceptional fit and reliable protection

This non-sterile glove features an advanced stretch formulation providing you with the ultimate comfort and fit. With its high strength and excellent protection features, this glove is ideal for nearly all clinical settings.

### KEY FEATURES AND BENEFITS:

- Advanced formulation for exceptional comfort and fit providing better tactile sensitivity
- An AQL better than industry standards delivers consistent, high quality production
- 200-count packaging reduces waste
- Available in size XXL
- Tested for use with over 55 chemicals and approved for use with chemotherapy drugs\*

\* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations. Users should test the suitability of this product against specific chemicals and the environment where used. Data for specific chemicals is available upon request.

**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used."

### Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Standard, moderate risk examination procedures
- Environmentally friendly facilities

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> Nitrile

Essential Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.07 / 2.8
Finger Thickness (mm/mil)	0.09 / 3.5	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), XXL (10.5 - 11)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Essential Series

## ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)	XXL (10.5 - 11)
Ansell Product Code	6034303	6034302	6034301	6034304	6034300	6034305

## PACKAGING AND STORAGE

Packaging	200 gloves per dispenser; 10 dispensers per case; 2000 gloves per case; *180 gloves per dispenser (Size XXL) * 1800 gloves per case (Size XXL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, ISO 9001

## Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com