



## MICRO-TOUCH® NITRA-TEX® EP™

### ALLERGY PREVENTION COMBINED WITH OUTSTANDING BARRIER PROTECTION AND CHEMICAL RESISTANCE.

Today, healthcare workers are highly aware of the risk of Type I allergies resulting from exposure to latex.

Made of 100% nitrile, Ansell's MICRO-TOUCH® Nitra-Tex® EP™ gloves are latex- and powder-free. Moreover, the 100% thiuram-free formulation and absence of chemicals significantly reduce the risk of Type IV allergies.

Nitrile gloves are stronger than vinyl gloves and most latex gloves. MICRO-TOUCH® Nitra-Tex® EP™ gloves have exceptionally good physical properties, as well as a strong resistance level to a large number of chemical products, including cytostatics.

Chlorination of the inner coating of the glove makes donning easy, while the absence of powder considerably reduces the risk of skin abrasion. Micro-textured finger surfaces ensure a good grip on instruments, and the gloves are also suitable for contact with foodstuffs.

The combination of easy donning, high comfort and sure grip make MICRO-TOUCH® Nitra-Tex® and MICRO-TOUCH® Nitra-Tex® EP™ the optimal choice of many healthcare professionals.

At 290 mm, MICRO-TOUCH® Nitra-Tex® EP™ offers additional length for extra protection.

This examination glove, CE Marked as Class I Medical Device, is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.

### PRODUCT DESCRIPTION

Material	Nitrile
Colour	Pale Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Textured
Internal Surface	Chlorinated

### SIZE / PRODUCT CODE

<b>x-small</b>	700121	<b>small</b>	700122
<b>medium</b>	700123	<b>large</b>	700124
<b>x-large</b>	700125		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

### Nitrile Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III  
 Registering Authority: Centexbel (0493): PPE  
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at [www.ansell.eu](http://www.ansell.eu).

#### Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J  
 Boulevard International 55, 1070 Brussels, Belgium  
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03  
<http://www.ansell.eu> • E-mail [info@ansell.eu](mailto:info@ansell.eu)

