

MICRO-TOUCH®

Examination Gloves

TECHNICAL DATASHEET

MICRO-TOUCH® Nitra-Tex® EP™

Nitrile Powder-Free Non-Sterile Examination Glove

Issue date 25-02-2020

PRODUCT DESCRIPTION

Material	Nitrile
Colour	Pale Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Textured
Internal Surface	Chlorinated

PHYSICAL PROPERTIES

Thickness (single)	Finger	0,155
Typical average values (mm)	Palm	0,120
	Cuff	0,100
Minimum length (mm)	290	
Strength (average values)	before ageing	after ageing
	595	595
Elongation at Break (%)	11,8	11,4
Force at break (N)	11,8	11,4

PACKAGING AND STORAGE

Packaging	10 x 100 pieces / 1000 pieces per case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight. Store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCT REFERENCES

Size / product code	x-small	700121	small	700122
	medium	700123	large	700124
	x-large	700125		


KEY FEATURES AND ADVANTAGES

- Additional length for extra protection
- Allergy prevention
- Comfort and ease-of-use
- Strong safety and chemical resistance


MANUFACTURING AND SAFETY STANDARDS

AQL (pinholes)	<ul style="list-style-type: none">• Waterleak tested• Conform European Norm EN 455-1: Inspection Level I AQL 1.5
Protein content	Not applicable: contains no natural rubber latex
Pyrogenicity	Absent of pyrogens as per ISO 10993-11
Primary skin irritation	Considered a negligible irritant as per ISO 10993-10
Skin Sensitization	No evidence of delayed dermal contact sensitisation as per ISO 10993-10
Viral Penetration	<ul style="list-style-type: none">• Passes Ansell's AVPP dynamic testing• Passes ASTM F 1671 using Phi X 174• Passes ISO 16604 using Phi X 174
Cytostatic Permeation	<ul style="list-style-type: none">• ACPPT tested under dynamic conditions• EN374-3 Breakthrough Times available
Sterilization	Not Applicable
CE Mark	Personal Protective Equipment: Category III Medical Device: Class I
Product Standards Conformance	EN 455 parts 1, 2, 3, 4 EN 420 EN ISO 374 part 1 EN 374 part 2 EN 16523 part 1 EN 374 part 4 EN ISO 374 part 5 EN 421
Manufacturing Standards	ISO 13485 ISO 9001
Food Contact	Regulation 1935/2004 Suitable for food contact EC
Registering Authority	Centexbel (0493): PPE

PICTOGRAMS



Please consult the Instruction for Use.



Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu