

# MICRO-TOUCH<sup>®</sup> NITRA-TEX<sup>™</sup>

Performance Series

## Nitrile, powder-free examination glove with strong barrier protection (only available in APAC)

This non-sterile glove is made from a softer synthetic formulation to offer a comfortable fit and feel. It also has an excellent resistance from a broad range of commonly used chemicals. High tear resistance ensures optimal protection, and a chlorinated inner coating facilitates ease of donning.

### KEY FEATURES AND BENEFITS:

- Excellent resistance to a broad range of commonly used chemicals
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs\*

\* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations


**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

\*\* This product is only available in Asia Pacific (APAC).

### Recommended For

- Tested for use with chemotherapy drugs
- Standard, moderate risk examination procedures
- Environmentally friendly facilities
- Protection from Type I latex allergy in HCW's or patients

### FEATURED TECHNOLOGIES



### MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> NITRA-TEX<sup>™</sup>

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.12 / 4.7
Finger Thickness (mm/mil)	0.15 / 5.9	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	250 / 9.8	Product Segmentation	Performance Series

## PACKAGING AND STORAGE

Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, EN 420:2003 + A1:2009, ISO 10993-10, ISO 13485, ISO 9001

## Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com