



MICRO-TOUCH® NITRILE ACCELERATOR-FREE

A UNIQUE COMBINATION OF BARRIER PROTECTION, ALLERGY PREVENTION AND COMFORT

Today, healthcare workers are highly aware of the risk of Type I allergies resulting from exposure to latex.

Made of 100% nitrile, Ansell's MICRO-TOUCH® Nitrile Accelerator-Free gloves contain no latex and are powder-free. Moreover, the complete absence of vulcanization accelerators including a 100% thiuram-free formulation significantly reduce the risk of Type IV allergies.

MICRO-TOUCH® Nitrile Accelerator-Free gloves provide effective barrier protection against punctures and a wide range of chemicals.

A thin glove, MICRO-TOUCH® Nitrile Accelerator-Free provides latex-like comfort. 'Thermo-elastic' material that adapts to the shape of the hand ensures even more comfort.

Chlorination of the inner surface of the glove facilitates donning, while the lack of powder considerably reduces the risk of skin abrasion. Micro-textured finger surfaces ensure a good grip on instruments. And the indigo blue colour makes easy product identification possible.

The combination of easy donning, high comfort and sure grip make MICRO-TOUCH® Nitrile Accelerator-Free the optimal choice of many healthcare professionals.

PRODUCT DESCRIPTION

| | |
|------------------|--------------|
| Material | Nitrile |
| Colour | Indigo |
| Shape | Ambidextrous |
| Cuff | Beaded |
| External Surface | Textured |
| Internal Surface | Chlorinated |

SIZE / PRODUCT CODE

| | | | |
|----------------|--------|--------------|--------|
| x-small | 700101 | small | 700102 |
| medium | 700103 | large | 700104 |
| x-large | 700105 | | |

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

Nitrile Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III
 Registering Authority: Centexbel (0493): PPE
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
 Boulevard International 55, 1070 Brussels, Belgium
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu