

# MICRO-TOUCH<sup>®</sup> Nitrile Accelerator-Free

Performance Series

## Nitrile, powder-free examination glove with advanced allergy prevention

This non-sterile accelerator-free offers a solution for those who suffer from chemical Type IV allergies. A comfortable fit and feel with micro-textured fingers ensure a secure and safe grip in wet or dry conditions.

### KEY FEATURES AND BENEFITS:


- Provides a high-level of barrier protection against punctures as well as excellent resistance from a broad range of commonly used chemicals.
- Chlorinated inner surface for easy donning
- Protects from latex Type I allergy
- 100% thiuram-free formulation significantly reduces the risk of Type IV allergies.

**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

### Recommended For

- Tested for use with chemotherapy drugs
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Sustained, high-risk examination procedures
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

### FEATURED TECHNOLOGIES



### MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> Nitrile Accelerator-Free

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.13 / 5.1	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Performance Series

## ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	700104	700103	700102	700105	700101

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

Category III, EN 374:2003, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 9001, ISO 10993-10:2002

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com