



MICRO-TOUCH® BLUE NITRILE

DON'T COMPROMISE ON RELIABILITY AND COMFORT

MICRO-TOUCH® Blue Nitrile examination gloves feature an advanced stretch formulation providing ultimate comfort and fit. With its high strength and excellent protection features, this glove is ideal for use in nearly all clinical settings.

MICRO-TOUCH® Blue Nitrile is a non-sterile, powder-free glove made of 100% nitrile. This examination glove does not contain natural rubber latex, thus addressing Type I (latex) allergy concerns. Nitrile gloves are recommended for conditions where high protection and chemical protection are required. This glove is tested for use with the following chemotherapy drugs: Cisplatin, Cyclophosphamide (Cytoxan), Dacarbazine (DTIC), Doxorubicin Hydrochloride, Etoposide (Toposar), Fluorouracil, Methotrexate, Mitomycin C, Paclitaxel (Taxol), Vincristine Sulfate.

This nitrile glove offers optimal 'touch and feel' and great fit to the user, instilling confidence to the user. Thanks to its nitrile formulation MICRO-TOUCH® Blue Nitrile has high elasticity and good memory allowing the nitrile film to adapt to the wearer's hand. The textured fingertips provide extra protection for wet or dry grip, and promote sure handling of materials and instruments.

PRODUCT DESCRIPTION

Material	Nitrile
Colour	Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Chlorinated

SIZE / PRODUCT CODE

x-small	313041060	small	313041065
medium	313041070	large	313041075
x-large	313041080		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

Nitrile Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III
 Registering Authority: Centexbel (0493): PPE

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
 Boulevard International 55, 1070 Brussels, Belgium
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu