



MICRO-TOUCH® ULTRA II

SUPERIOR COMFORT AND SAFETY

Like a second skin, Ansell's MICRO-TOUCH® Ultra II glove offers comfortable, precise handling with excellent sensitivity.

The inner surface of MICRO-TOUCH® Ultra II gloves is chlorinated to ensure easy donning, and it has a smooth external surface that makes it ideal for patient handling and facilitates double gloving.

Moreover, the powder-free design and low protein content minimise the risk of latex-related allergies for the wearer: extensive leaching phases and post-production washing cycles ensure MICRO-TOUCH® Ultra II has very low residual protein content. This (and the fact that gloves are powder-free) significantly reduces the risk of Type I allergies. The glove also contains a low level of residual process chemicals, minimising the risk of Type IV allergies.

A reference in the powder-free medical environment, MICRO-TOUCH® Ultra II offers a secure barrier with a high degree of tear resistance and is the practical choice, combining tactility, ease of donning and protection.

This examination glove, CE Marked as Class I Medical Device, is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.

PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Beige
Shape	Ambidextrous
Cuff	Beaded
External Surface	Textured, Chlorinated
Internal Surface	Chlorinated

SIZE / PRODUCT CODE

x-small	553691	small	553692
medium	553693	large	553694
x-large	553695		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

Natural Latex Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III
 Registering Authority: Centexbel (0493): PPE
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu