

MICRO-TOUCH[®] Ultra II

E s s e n t i a l S e r i e s

Latex, powder-free, examination glove optimal for tactility, ease of donning and protection

This non-sterile glove has a chlorinated inner surface designed to ensure easy and fast donning. Its fully textured surface provides excellent tactile sensitivity, making it ideal for patient handling and double gloving. Ansell's natural rubber latex formulation provides a comfortable, well-fitting glove while extensive leaching phases and post-production washing cycles ensure a very low residual protein content. CE Marked as a Class I Medical Device, it is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.


KEY FEATURES AND BENEFITS:

- Unique latex formulation offers comfortable, precise handling with excellent sensitivity
- Chlorinated inner coating to ensure easy donning
- Secure barrier protection with a high degree of tear resistance

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Brief, low-risk examination procedures



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Ultra II

Essential Series

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.14 / 5.5
Finger Thickness (mm/mil)	0.16 / 6.3	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5
Product Segmentation	Essential Series		

ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	553694	553693	553692	553695	553691

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Category III, EN 374:2003, EN 374:2003, EN 420:2003 + A1:2009, EN 421:2010, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10993-10:2002, ISO 13485, ISO 9001

Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)	Zinc Diethyldithiocarbamate (ZDEC)
Zinc Mercaptobenzothiozole (ZMBT)	

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com