



PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Beige
Shape	Ambidextrous
Cuff	Beaded
External Surface	Textured
Internal Surface	Polyacrylate coated

MICRO-TOUCH® STERILE

SUPERIOR COMFORT AND GRIP WHEN STERILITY IS REQUIRED

For procedures requiring sterility and when surgical gloves are not required, Ansell has developed MICRO-TOUCH® Sterile examination gloves.

To ensure sterility is preserved, they are packaged in pairs with a protective double wrapping, clearly indicating the glove size.

Users find MICRO-TOUCH® Sterile a very comfortable glove, and donning is easy even with damp hands. An interior polyacrylate coating makes the glove even softer, while the textured fingers provide a safe grip when handling instruments, even with damp hands.

Powder-free and with low levels of protein and residual chemicals, they minimise the risk of Type I allergies.

Ansell MICRO-TOUCH® Sterile gloves are the sensible solution for sterile examination (non-surgical) procedures.

SIZE / PRODUCT CODE

small	553342	medium	553343
large	553344		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

Natural Latex Powder-Free Sterile Examination Glove

MD Class I sterile
If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. Sterile unless this package is opened or damaged. This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
Boulevard International 55, 1070 Brussels, Belgium
Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu

