



## MICRO-TOUCH® HYDRACARE®

### A UNIQUE COMBINATION OF PROTECTION AND SKIN CARE

To minimise the effects of continuous glove use on hands, Ansell offers MICRO-TOUCH® HydraCare®.

This powder-free latex examination glove contains an interior glycerol moisturiser that helps maintain skin suppleness. This has a proven corrective effect on defective skin, restoring hydration and elasticity.

Single chlorination combined with silicone and the moisturiser guarantee easy donning, even with damp hands.

A pale green colour ensures easy product identification, while the glove's pleasant scent is much appreciated by healthcare workers and patients alike.

Extensive leaching phases and post-production washing cycles ensure the glove has a very low residual protein content. This, along with the fact that it is powder-free, significantly reduces the risk of Type I allergies. Furthermore, the low level of residual process chemicals and absence of thiurams minimise the risk of Type IV allergies.

MICRO-TOUCH® HydraCare® offers high-quality protection, comfort and ease of use.

### PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Pale Green
Shape	Ambidextrous
Cuff	Beaded
External Surface	Micro-Textured
Internal Surface	Slightly chlorinated, Treated with silicone and moisturizer

### SIZE / PRODUCT CODE

<b>x-small</b>	3030120XS	<b>small</b>	30301200S
<b>medium</b>	30301200M	<b>large</b>	30301200L
<b>x-large</b>	3030120XL		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

### Natural Latex Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III  
 Registering Authority: Centexbel (0493): PPE  
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at [www.ansell.eu](http://www.ansell.eu).

#### Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J  
 Boulevard International 55, 1070 Brussels, Belgium  
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03  
<http://www.ansell.eu> • E-mail [info@ansell.eu](mailto:info@ansell.eu)

