

# MICRO-TOUCH<sup>®</sup> HYDRACARE<sup>™</sup>

Performance Series

## Latex, powder-free examination glove with a skin moisturizing coating

This non-sterile glove contains an interior glycerol moisturizer that helps maintain skin suppleness, restoring hydration and elasticity. In combination with the inner chlorinated surface, this glove guarantees easy donning, even with damp hands. Ansell's natural rubber latex formulation provides a comfortable, well-fitting glove while extensive leaching phases ensure a very low residual protein content. This glove is CE Marked as a Class I Medical Device, ensuring more consistent quality and providing protection for a wide range of applications.

### KEY FEATURES AND BENEFITS:

- Fully textured surface provides extra protection for wet or dry grip and against slippage
- Glycerol inner coating actively retains moisture and rehydrates skin
- Excellent elasticity and memory for a secure, comfortable fit

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

### Recommended For

- Enhanced grip in wet environments
- Standard, moderate risk examination procedures
- Extra comfort during long, sustained procedures
- Damp donning
- Helping protect skin from drying or cracking

### FEATURED TECHNOLOGIES

**HYDRACARE<sup>™</sup>**  
Inner Coating Technology



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> HYDRACARE<sup>™</sup>

Performance Series

## PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	Green	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.13 / 5.1
Finger Thickness (mm/mil)	0.15 / 5.9	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5
Product Segmentation	Performance Series		

## ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	553334	553333	553332	553335	553331

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *90 gloves per dispenser (Size XL) *900 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D3578, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10993-10:2002, ISO 13485, ISO 10993-5, ISO 9001, ISO 10993-10:2003

## Sensitisers

Zinc Diethyldithiocarbamate (ZDEC)	Zinc Dibutyldithiocarbamate (ZDBC)
Zinc Mercaptobenzothiozole (ZMBT)	

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com