



## MICRO-TOUCH® COATED

### EASY DONNING AND MAXIMUM COMFORT

Thanks to a polyacrylate internal coating which ensures easy and fast donning, Ansell's MICRO-TOUCH® Coated gloves are ideal for users who must switch gloves quickly and frequently.

Gloves are very comfortable and a range of five sizes ensures a perfect fit.

This examination glove, CE Marked as Class I Medical Device, is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.

The latex formulation and optimised film thickness offer high levels of elasticity, flexibility and sensitivity.

Extensive leaching phases and post-production washing cycles ensure the glove has a very low residual protein content. This, along with the fact that it is powder-free, significantly reduces the risk of Type I allergies. Furthermore, the low level of residual process chemicals and absence of thiurams minimise the risk of Type IV allergies.

With a low allergenic profile, MICRO-TOUCH® Coated is the practical choice for healthcare professionals who want to establish or maintain a powder-free environment.

### PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Beige
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Polyacrylate coated

### SIZE / PRODUCT CODE

<b>x-small</b>	553301	<b>small</b>	553302
<b>medium</b>	553303	<b>large</b>	553304
<b>x-large</b>	553305		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

### Natural Latex Powder-Free Non-Sterile Examination Glove

MD Class I & PPE Cat. III  
 Registering Authority: Centexbel (0493): PPE  
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at [www.ansell.eu](http://www.ansell.eu).

#### Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J  
 Boulevard International 55, 1070 Brussels, Belgium  
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03  
<http://www.ansell.eu> • E-mail [info@ansell.eu](mailto:info@ansell.eu)

