

MICRO-TOUCH®

Examination Gloves

TECHNICAL DATASHEET

MICRO-TOUCH® Coated

Natural Latex Powder-Free Non-Sterile Examination Glove

Issue date 25-02-2020

PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Beige
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Polyacrylate coated

PHYSICAL PROPERTIES

Thickness (single)	Finger	0,150
Typical average values (mm)	Palm	0,110
	Cuff	0,080
Minimum length (mm)	245	
Strength (average values)	before ageing	after ageing
	715	620
Elongation at Break (%)	10,0	8,9

PACKAGING AND STORAGE

Packaging	10 x 100 / 1000 per case, 10 x 80 / 800 per case (size:x-large)
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight. Store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCT REFERENCES

Size / product code	x-small 553301	small 553302
	medium 553303	large 553304
	x-large 553305	


KEY FEATURES AND ADVANTAGES

- Easy donning
- Excellent comfort and sensitivity
- High level of safety and low allergenicity

MANUFACTURING AND SAFETY STANDARDS

AQL (pinholes)	<ul style="list-style-type: none">• Waterleak tested• Conform European Norm EN 455-1: Inspection Level I AQL 1.5
Protein content	Very low level: Typical results < 30 µg/g for Modified Lowry & HLPC (EN 455-3)
Pyrogenicity	Absent of pyrogens as per ISO 10993-11
Primary skin irritation	Considered a negligible irritant as per ISO 10993-10
Skin Sensitization	No evidence of delayed dermal contact sensitization as per ISO 10993-10
Viral Penetration	<ul style="list-style-type: none">• Passes ASTM F 1671 using Phi X 174• Passes ISO 16604 using Phi X 174
Sterilization	Not Applicable
CE Mark	Personal Protective Equipment: Category III Medical Device: Class I
Product Standards Conformance	EN 455 parts 1, 2, 3, 4 EN 420 EN ISO 374 part 1 EN 374 part 2 EN 16523 part 1 EN 374 part 4 EN ISO 374 part 5 EN 421
Manufacturing Standards	ISO 13485 ISO 9001
Registering Authority	Centexbel (0493): PPE

PICTOGRAMS



Please consult the Instruction for Use.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu