

MICRO-TOUCH[®] Coated

Essential Series

Latex, powder-free, examination glove coated for quick and frequent donning

This non-sterile glove features a polyacrylate inner coating ensuring easy and fast donning, making it ideal for users who must switch gloves frequently. Ansell's natural rubber latex formulation provides a comfortable, well-fitting glove while extensive leaching phases ensure a very low residual protein content. CE Marked as a Class I Medical Device, it is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications.

KEY FEATURES AND BENEFITS:

- Polyacrylate internal coating facilitates superior ease of donning
- Latex formulation and optimized film thickness offer high levels of elasticity, flexibility and sensitivity
- Low allergenicity latex reduces risk of contact dermatitis due to absence of thiurams

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

Recommended For

- Brief, low-risk examination procedures
- Damp donning



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Coated

Essential Series

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.11 / 4.3
Finger Thickness (mm/mil)	0.14 / 5.5	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5
Product Segmentation	Essential Series		

ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	553304	553303	553302	553305	553301

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Category III, EN 374:2003, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 421:2010, FDA21 CFR 177-2600, ISO 11193-1:2008, ISO 13485, ISO 9001, ISO 10993-10:2010

Sensitisers

Zinc Diethyldithiocarbamate (ZDEC)	Zinc Dibutyldithiocarbamate (ZDBC)
Zinc Mercaptobenzothiozole (ZMBT)	

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com