



PRODUCT DESCRIPTION

Material	Nitrile
Colour	Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Chlorinated

MICRO-TOUCH® NITRILE STERILE

NITRILE, STERILE PAIR PACKED, POWDER-FREE EXAMINATION GLOVE FEATURING EXTENDED CUFF

This sterile glove is made from a synthetic, nitrile formulation that provides high elasticity with superior barrier protection and puncture resistance. It is an ideal solution for a range of sterile examination (non-surgical) procedures.

KEY FEATURES & BENEFITS

- Packed as a pair of gloves per pack
- 12" extended cuff for additional protection
- Textured fingers provide extra grip and reliability
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs*

* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

RECOMMENDED FOR:

- Protection from latex Type I allergy in HCW's or patients
- Sterile procedures outside of the O.R.
- Tested for use with chemotherapy drugs

SIZE / PRODUCT CODE

small	6034151	medium	6034152
large	6034153	x-large	6034154

Nitrile Powder-Free Sterile Examination Glove

MD Class I & PPE Cat. III
 Registering Authority: Centexbel (0493): PPE
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. Sterile unless this package is opened or damaged.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
 Boulevard International 55, 1070 Brussels, Belgium
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu