



### **Product description**

Order number	72.694.006
Product description	Screw Cap Micro Tube, 2 ml, PP, with skirted base, with knurls, with assembled cap, with printed writing space and graduation, sterile, 100 pcs./bag

### **Product characteristics**

Label/ Print	With print
Colour of print/label	white
Graduation	yes
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Diameter	10,8 mm
Length of product	44 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### Purity & certification

Satisfies the requirement	IATA
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Purity standard	sterile
Sterilisation	yes
Pyrogen-free	no
Lot or serial number	Lot number

### Packaging

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0,0135 cbm
Weight of case	1,85 kg
EAN of inner box	4038917117180
EAN of outer box	4038917015486

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

