



### **Product description**

Order number	75.562.011
Product description	Container with screw cap 100 ml, PP brown, with enclosed assembled brown cap

### **Product characteristics**

Graduation	yes
Type of cap	Screw cap
Light protection	yes
Application	Collection and storage of urine

### **Size**

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	brown

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	250
Pcs. / inner box	50
Pcs. / box	250
Pcs. / pallet	7500
Depth of case	443 mm
Width of case	323 mm
Height of case	330 mm
Case volume	0,0472 cbm
Weight of case	4,59 kg
EAN of inner box	4038917113618
EAN of outer box	4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

