




UTM™: UNIVERSAL
TRANSPORT MEDIUM
OFTEN IMITATED,
NEVER REPLICATED.

COLLECTION | TRANSPORT | PRESERVATION



> UTM™ ADVANTAGE

- ROOM TEMPERATURE STABLE
- RAPID RELEASE AND DISPERSION of patient sample material and virus particles during vortexing due to three glass beads in each tube.
- UNIQUE MEDIA FORMULATION includes antibiotics to inhibit bacteria and fungal flora in patient samples
- SAFE AND RELIABLE MOLDED BREAKPOINT
- SAFE, SHATTERPROOF, STAND UP tubes with distinctive internal conical shape enabling centrifugation of samples.
- AVAILABLE WITH REGULAR AND MINITIP FLOQSWABS™
- COST EFFICIENT large 3mL media fill volume allows for multiple tests on the same specimen. 1mL fill volume prevents sample dilution.

For collection, transport,
maintenance and long term freeze
storage of Virus, Chlamydia,
Mycoplasma and Ureaplasma

Available as filled tubes in bulk or as patient sample collection packs which can incorporate various combinations of filled tubes with either FLOQSwabs™ (regular or mini tip) or polyester fiber swabs.



UTM™ COPAN UTM™ ORDERING INFORMATION

KITS OF UTM™ WITH FLOCKED SWABS



Cat. No.	Tube Size	Swab Type	Product Description	Unit of Measure (UOM)
306C	16x100mm	Flocked	Conical tube filled with 3mL UTM™ medium packaged with one regular FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box
307C	16x100mm	Flocked	Conical tube filled with 3mL UTM™ medium packaged with one minitip FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box
305C	16x100mm	Flocked	Conical tube filled with 3mL UTM™ medium packaged with one flexible minitip FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box
359C	12x80mm	Flocked	Conical tube filled with 1mL UTM™ medium packaged with one regular FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box
361C	12x80mm	Flocked	Conical tube filled with 1mL UTM™ medium packaged with one minitip FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box
360C	12x80mm	Flocked	Conical tube filled with 1mL UTM™ medium packaged with one flexible minitip FLOQSwab™, sterile	50 kits per shelf pack 6 x 50 kits per box

KITS OF UTM™ WITH POLYESTER SWABS



Cat. No.	Tube Size	Swab Type	Product Description	Unit of Measure (UOM)
302C	16x100mm	Polyester	Conical tube filled with 3mL UTM™ medium packaged with 2 regular size swabs with plastic applicator and polyester tip	50 kits per shelf pack 6 x 50 kits per box
340C	16x100mm	Polyester	Conical tube filled with 3mL UTM™ medium packaged with 1 regular size swab with plastic applicator and polyester tip and 1 minitip size swab with wire applicator and polyester tip	50 kits per shelf pack 6 x 50 kits per box

Additional packaging configurations may be available. Contact COPAN or your authorized COPAN distributor.

UTM™ TUBES IN BULK



Cat. No.	Tube Size	Product Description	Unit of Measure (UOM)
330C	16x100mm	Conical tube filled with 3mL UTM™ medium	50 tubes per shelf pack 6 x 50 tubes per box
331C	25x90mm	Conical tube filled with 10mL UTM™ medium	15 tubes per shelf pack 6 x 15 tubes per box
350C	12x80mm	Conical tube filled with 1mL UTM™ medium	50 tubes per shelf pack 6 x 50 tubes per box

SWAB KITS FOR USE WITH UTM™

Cat. No.	Product Description	Unit of Measure (UOM)
340CS01	One regular size swab with plastic applicator and polyester tip and one minitip size swab with wire applicator and polyester tip	100 kits per shelf pack 10 x 100 kits per box
302CS01	Two regular size swabs with plastic applicator and polyester tip	100 kits per shelf pack 10 x 100 kits per box

COPAN DIAGNOSTICS, INC.

26055 JEFFERSON AVENUE
MURRIETA, CA 92562 USA

TL. FREE 800.216.4016
TEL 951.696.6957
FAX 951.600.1832



Innovating Together™

UTM09151C