

GAMMEX® Latex Powder-Free with AMT™

S U R G I C A L G L O V E S

Latex, powder-free surgical glove with a proprietary antimicrobial inner coating for active barrier protection

Realised through extensive research and development, Ansell's proprietary AMT™ Technology enables the application of chlorhexidine gluconate (CHG) to the glove, providing rapid anti-viral and anti-bacterial action. This glove suppresses the regrowth of bacteria on healthcare workers' hands and delivers protection even when a glove breach goes unnoticed, decreasing the risk of surgical site infection.¹

KEY FEATURES & BENEFITS

- Proprietary AMT Technology™ enables the application of chlorhexidine gluconate (CHG) to the inner coating, suppressing the users' regrowth of bacteria while providing rapid anti-viral and anti-bacterial action for added protection even when a glove breach occurs
- Kills at least 99% of Hep-C virus and HIV-1 strain Mn as early as 1 minute²
- Kills between 99 and 99.999% of eight common gram-positive, gram-negative and drug-resistant bacteria as early as 1 minute²
- This antimicrobial glove contains moisturizers, emollients and anti-irritant components ensuring a skin friendly glove while providing enhanced protection

NOTE: Not proven to protect against blood-borne infections where the skin is broken, cut or punctured.

1. Suppression of surgeons' bacterial hand flora during surgical procedures with a new antimicrobial surgical glove. Assadian O, Kramer A, Ouriel K, Suchomel M, McLaws ML, Rottman M, Leaper D, Assadian A. Surg Infect. 2014 Feb;15(1):43-9.

2. Data on file. HIV-1 strain Mn is a well-studied, virulent strain of HIV-1.

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

RECOMMENDED FOR

- Damp donning
- Extra comfort during long, sustained procedures
- Use as an outer glove or underglove when double gloving

FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex Powder-Free with AMT[™]

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural rubber latex
Powder Content	Powder-free
Color	Natural
Shape	Anatomic with straight fingers
Internal Glove Surface	Antimicrobial coated
External Glove Surface	Smooth with micro-textured finish
Cuff Style	Straight with SUREFIT [™] Technology
Thickness	Standard
Grip Level	Moderate
Allergy Prevention	None
Double Gloving Recommendation	Outer glove or underglove
Freedom From Holes (AQL)	In process control before packing: 0.4 (Inspection level GI), Manufacturing final release: 1.0 (Inspection level GI)
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978

ORDERING INFORMATION

Size	5 ½	6	6 ½	7	7 ½	8	8 ½	9
Ansell Product Code	353441	353442	353443	353444	353445	353446	353447	353448

PACKAGING AND STORAGE

Packaging	25 pairs per box; 4 boxes per carton/case; 100 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Product Standards	EN 455 Part 1-4, EN ISO 374-1, EN 374-2 and -4, EN 16523-1, EN ISO 374-5, EN 420, ISO 10282, AS/NZS 4179, JIS T9107
Quality / Environmental Standards	ISO 9001, ISO 13485, ISO 14001, ISO 11137-Part 1, EN 556

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com