



PurFlock Ultra® Flocked Swab - Dry Transport System

Ultrafine tip, polystyrene handle with breakpoint at 100 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.14" Dia. (3.556 mm)
- 0.66" Length (16.764 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.036" (0.914 mm)
- 5.99" Length (147.62 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Rapid diagnostic testing
- Forensic evidence collection
- Specimen collection

Ref #	Packaging
25-3317-U BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®