

C13600 **Viscose Universal Cap**
C13700 **Medical Device Class I**

FEATURES & BENEFITS

- Very extensible aerated nonwoven
- Very aesthetic and covering cap
- Way of use:
 1. To be aesthetic, seams can be returned.
 2. Adjust the cap by rolling up nonwoven from the bottom.



Recommendation:

- Can be used by nurse, surgeon or for patient in operating room

TECHNICAL DESCRIPTION

Product name: 100 % Viscose universal cap
 Product type: Single use, non sterile
 Material : 80 % Viscose + 20 % Polymer binder
 Confection : Sewn
 Color : Blue / Green
 Dimensions : Circumference : 56 cm (± 2 cm)
 Height : 32 cm (± 2 cm)
 Poids : 7,2 g (± 10%)
 Packaging : Dispenser box of 100 units
 Origin : France



| Item | Color |
|--------|---|
| C13600 | Blue  |
| C13700 | Green  |



Mechanical information (Measurements done by KOLMI Laboratory)

| TEST | Standard | | Average | Unit | Tolerance |
|------------------|------------------|---------|---------|------------------|-----------|
| Grammage | NF EN ISO 9073-1 | | 40 | g/m ² | ± 10% |
| Thickness | NF EN ISO 9073-2 | | 1,2 | mm | ± 10% |
| Tensile strength | NF EN ISO 9073-3 | Dry MD: | 35 | N | ± 10% |
| | | Dry CD: | 4 | N | ± 10% |

Manufacturing Microbiological Information (by MICROSEPT laboratory)

Bioburden evaluation meets the requirements of NF EN 11737:2006 Standard
 Additional microbiological controls: ASR, E.coli, Staphylococcus.

C13600
C13700

Viscose Universal Cap
Medical Device Class I

STANDARDS & CERTIFICATIONS

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC
Manufactured under an ISO 13485 certified system
Complies with the harmonized standard applicable to EN ISO 14971.

STORAGE CONDITIONS

Normal conditions of conversation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

LOGISTICAL INFORMATION



Outer case
specifications

| Item | Size mm | Gross weight | QTY/ pallet |
|------------------|-----------------|-----------------|-----------------------------------|
| C13600 C13700 | 575 x 225 x 370 | 4,6 | 35 cases (5 layers of 7 cases) |

Dispenser box
specifications

| Size mm | QTY |
|-----------------|-----------------------------------|
| 200 x 108 x 350 | 5 dispenser boxes of 100 units |

 EAN Code

| Item | Case | Dispenser box |
|--------|-------------------|-------------------|
| C13600 | 3 662 036 001 108 | 3 662 036 000 088 |
| C13700 | 3 662 036 001 115 | 3 662 036 000 095 |

