

ENCORE® Non-Latex PI Underglove

S U R G I C A L G L O V E S

Polyisoprene, skin-friendly surgical underglove, an ideal non-latex double gloving solution

KEY FEATURES AND BENEFITS:

- Ultra-thin underglove offering optimal sensitivity and comfort when double gloving
- Green color for easy glove breach identification when double gloving
- Softer and 'stretchier' PI glove featuring superior comfort compared to Ansell's latex gloves
- Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities for optimized allergenic profile
- Smooth finish offers adequate grip, allowing the glove to be used in a single glove environment

**Data on file: HRIPT (Human Repeat Insult Patch Testing) conducted by an independent laboratory in the U.S.*

PI products featuring PI-KARE Technology may not yet be available in your country. Please contact your Ansell representative.

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Use as an underglove when double gloving
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

FEATURED TECHNOLOGIES



ENCORE® The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



ENCORE[®] Non-Latex PI Underglove

SURGICAL GLOVES

PRODUCT INFORMATION

| | |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------|
| Material | PI (Polyisoprene) |
| Powder Content | Powder-Free |
| Color | Green |
| Thickness | Micro |
| Cuff Style | Beaded with SUREFIT [™] Technology |
| External Glove Surface | Smooth with Micro-Textured Finish |
| Internal Glove Surface | Polymer Coated with DERMASHIELD [™] Technology |
| Allergy Prevention | Latex (Type I) |
| Grip Level | Low |
| Double Gloving Recommendation | Underglove |
| Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k) |
| Freedom from Holes (Inspection level I) | 0.65 AQL |
| Product Standards | ASTM D6978, EN 16523-1, EN 374-2 and -4, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016 |
| Quality / Environmental Standards | EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001 |
| Packaging | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case |
| Shelf Life | 3 years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Vulcanization Chemical Accelerators | Xanthates, Zinc Dialkyldithiocarbamate |

ORDERING INFORMATION

| Size | 5.5 | 6 | 7 | 7.5 | 8 | 8.5 | 9 |
|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Ansell Product Code | 330114055 | 330114060 | 330114070 | 330114075 | 330114080 | 330114085 | 330114090 |

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21)2407 2286/88/62
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2020 Ansell Limited. All Rights Reserved.