



ERGONOMIC SURGICAL GLOVES



ENCORE® NON-LATEX PI UNDERGLOVE

DELIVERING EXCEPTIONAL PERFORMANCE AND PROTECTION FOR EVERY HAND

Ansell is committed to delivering innovative products that address patient and healthcare worker safety in the complete perioperative environment - before, during, and after surgery.

Because of the significant impact of latex and chemical allergies, we have focused on developing innovative products that incorporate new non-latex materials and technologies to bring you the natural fit and comfort you love in ENCORE® gloves.

The ENCORE® Non-Latex product line offers a very complete range of non-latex options, taking allergy protection beyond latex with advanced protection against chemical allergic reactions for healthcare workers.

ENCORE® Non-Latex PI UnderGlove is a DPG-free polyisoprene option for healthcare workers suffering from latex allergies and chemical allergies caused by DPG.

A non-latex, DPG-free underglove offering optimal protection.

Its green colour makes it ideal as a non-latex double-gloving solution, together with the white ENCORE® Non-Latex PI Hybrid or ENCORE® Ultra.

PRODUCT DESCRIPTION

Material	Synthetic Polyisoprene
Colour	Green
Shape	Anatomic
Cuff	Beaded cuff with adhesive band
External Surface	Smooth, Chlorinated & Siliconized
Internal Surface	Polyurethane coated & Siliconized

SIZE / PRODUCT CODE

5.5	340114055	6	340114060
6.5	340114065	7	340114070
7.5	340114075	8	340114080
8.5	340114085	9	340114090

ENCORE® gloves are designed and engineered for surgeons, OR nurses and surgical technicians seeking protection and safety with advanced levels of ergonomics, dexterity and comfort during surgery. ENCORE® gloves afford the kind of hand protection that elevates users' confidence and provide optimum comfort during surgery.

Synthetic Polyisoprene Powder-Free Sterile Surgical Glove

MD Class IIa
Registering Authority: British Standards Institution (2797): MD
If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. Sterile unless this package is opened or damaged.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
Boulevard International 55, 1070 Brussels, Belgium
Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu

