



ERGONOMIC SURGICAL GLOVES



ENCORE® NON-LATEX®

OUTSTANDING PROTECTION FOR USERS WITH TYPE I (LATEX) AND TYPE IV (CHEMICAL ACCELERATOR) ALLERGIES

ENCORE® Non-Latex® is a powder-free surgical glove made of Neoprene.

Accelerator-free, the new SENSOPRENE® formulation delivers an exceptional surgical glove offering unprecedented comfort and sensitivity with advanced allergy protection.

No chemicals (thiurams, MBT, carbamates) known to contribute to Type IV allergies are used, making it safe for latex-sensitive (Type I) and chemical-sensitive (Type IV) healthcare professionals and patients.

New ENCORE® Non-Latex® gloves protect healthcare workers and patients against additional causes of skin reactions with features such as casein-free formulation - to minimise reactions to animal proteins - and CPC-free polymer coating – thus eliminating irritation caused by cetylpyridinium chloride (CPC).

With its new creamed colour, it offers enhanced tactile sensitivity in the fingertips and allows for easy and comfortable double gloving.

Offering Ansell's SureFit™ technology to prevent cuff roll-down, ENCORE® Non-Latex® provides healthcare professionals with superior barrier protection.

PRODUCT DESCRIPTION

Material	Neoprene
Colour	Cream
Shape	Anatomic
Cuff	SureFit™ technology: Beaded cuff with adhesive band
External Surface	Textured, Chlorinated & Siliconized
Internal Surface	Polyurethane coated & Siliconized

SIZE / PRODUCT CODE

5.5	340111055	6	340111060
6.5	340111065	7	340111070
7.5	340111075	8	340111080
8.5	340111085	9	340111090

ENCORE® gloves are designed and engineered for surgeons, OR nurses and surgical technicians seeking protection and safety with advanced levels of ergonomics, dexterity and comfort during surgery. ENCORE® gloves afford the kind of hand protection that elevates users' confidence and provide optimum comfort during surgery.

Neoprene Powder-Free Sterile Surgical Glove

MD Class IIa
Registering Authority: British Standards Institution (2797): MD, Centexbel (0493): PPE
If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. Sterile unless this package is opened or damaged.

Ref-XXXX-2020607

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
Boulevard International 55, 1070 Brussels, Belgium
Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu

