

schülke -†



The compact and practical solution for flexible disinfection columns.

schülke disinfection column

Our Plus:

- compatible as standard with hyclick®, sm 2 and KH dispensers
- perfect for entrance areas and places where wall mounting of dispensers is not desired or possible
- powder-coated steel for optimal and easy cleaning
- simple and quick installation

Description


The schülke disinfection column is a space-saving, stationary column for hand disinfectant dispensers 500 ml and 1000 ml.

Technical Data

Material	Steel, white powder-coated
Weight:	16,5 kg

Order information

Column	Art.-No.
schülke disinfection column	60000141



schülke disinfection column

schülke disinfection column with hyclick®, sm 2, KH dispensers



hyclick®




sm 2



KH

Drilling pattern for hyclick®, sm 2, KH, KHN touchless dispenser



- 1 - sm 2 500
- 2 - sm 2 universal + KH 500 + hyclick®
- 3 - KH 1000 + KHN touchless 500/1000
- 4 - drip tray SH 1

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes and in compliance with high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:
Customer Sales Service
Phone: +49 (0)40 521 00-666
E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt,

Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com

