





## schülke instrument-bathes

schülke offers compatible bathes for every requirement. In all sizes, including white or transparent lid, with instrument-box and drainage tap.



	Article	Article-code	Dimensions* inside (in mm)	Dimensions* outside (in mm)
	<b>3 L bath Complete - Set</b>			
	lid white	144 307	250 x 150 x 120	300 x 185 x 125
	lid transparent	144 308	250 x 150 x 120	300 x 185 x 125
	<b>5 L bathes Complete - Set</b>			
	lid white	144 407	500 x 170 x 120	555 x 205 x 135
	lid transparent	144 408	500 x 170 x 120	555 x 205 x 135
	<b>10 L bathes Complete - Set</b>			
	lid white	144 507	340 x 260 x 170	400 x 300 x 180
	lid transparent	144 508	340 x 260 x 170	400 x 300 x 180
	<b>30 L bathes Complete - Set</b>			
	lid white / including drainage tap	144 609	560 x 370 x 225	620 x 405 x 240
	lid transparent / including drainage tap	144 610	560 x 370 x 225	620 x 405 x 240
	lid white / without drainage tap	144 607	560 x 370 x 225	620 x 405 x 240
	lid transparent / without drainage tap	144 608	560 x 370 x 225	620 x 405 x 240

Several single components are available on request. Reprocessing not over 55 °C.

\* dimensions: length x width x height



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).

**Schülke & Mayr GmbH**  
22840 Norderstedt, Germany  
Phone +49 (0) 40 - 521 00 - 0  
Fax +49 (0) 40 - 521 00 - 318  
www.schuelke.com  
mail@schuelke.com

**Schülke & Mayr UK Ltd.**  
1 Jenkin Road  
GB-Sheffield S9 1AT  
Phone +44 (0) 1142 - 54 35 - 00  
Fax +44 (0) 1142 - 54 35 - 01  
mail.uk@schuelke.com

**Schülke & Mayr Ges.m.b.H**  
Seidengasse 9  
1070 Wien, Austria  
Phone +43 (0) 1 - 523 25 01 - 0  
Fax +43 (0) 1 - 523 25 01 - 60  
office.austria@schuelke.com

2361 / 1 / 03.10 / A / westwerk - This product information is not automatically updated.