



ENCORE NON-LATEX PI HYBRID IN SMART PACK

A NON-LATEX, DPG-FREE GLOVE PERFECT FOR GENERAL SURGERY

ENCORE® Non-Latex PI Hybrid is a powder-free surgical glove made from a unique blend of synthetic lattices (polyisoprene & polychloroprene). It is free from latex proteins making it safe for latex-sensitive (Type I) healthcare professionals and patients.

ENCORE® Non-Latex PI Hybrid is a DPG-free option for healthcare workers suffering from latex allergies and chemical allergies caused by DPG, an increasing problem.

This new glove protects healthcare workers and patients against additional causes of skin reactions with features such as casein free formulation - to minimize reactions to animal proteins - and CPC-free polymer coating - thus eliminating irritation caused by cetylpyridinium chloride (CPC).

Manufactured using a special former, ENCORE® Non-Latex PI Hybrid has a curved fit and a greater width across the hand: this reduces pressure on overworked areas such as the ball of the thumb and palm. The result is good comfort for the wearer, which is especially valuable during lengthy operative procedures. Ansell's SureFit™ technology helps prevent cuff roll-down for a secure & comfortable fit.

Its cream colour makes it ideal as a non-latex double-gloving solution, together with the green ENCORE® Non-Latex PI Underglove.

PRODUCT DESCRIPTION

Material	Synthetic Polyisoprene and Neoprene
Colour	Cream
Shape	Anatomic
Cuff	SureFit™ technology: Beaded cuff with adhesive band
External Surface	Textured, Chlorinated & Siliconized
Internal Surface	Polyurethane coated & Siliconized

SIZE / PRODUCT CODE

5.5	340105055	6	340105060
6.5	340105065	7	340105070
7.5	340105075	8	340105080
8.5	340105085	9	340105090

ENCORE® gloves are designed and engineered for surgeons, OR nurses and surgical technicians seeking protection and safety with advanced levels of ergonomics, dexterity and comfort during surgery. ENCORE® gloves afford the kind of hand protection that elevates users' confidence and provide optimum comfort during surgery.

Synthetic Polyisoprene and Neoprene Powder-Free Sterile Surgical Glove

MD Class Ila
 Registering Authority: British Standards Institution (2797): MD, Centexbel (0493): PPE
 If relevant, read carefully the instructions for use flyer accompanying the product. Single use only. Sterile unless this package is opened or damaged.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
 Boulevard International 55, 1070 Brussels, Belgium
 Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu