



### **Product description**

Order number	70.3040
Product description	Pipette tip, 300 µl, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Brand and Finnpiquette and products of identical design, 1000 pcs./bag

### **Product characteristics**

Graduation	yes
Volume of work	300 µl
Filter	no
Coding colour	orange
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	10000
Type of inner packaging	bag
Pcs. / inner box	1000
Pcs. / box	10000
Pcs. / pallet	420000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0,0396 cbm
Weight of case	4,18 kg
EAN of inner box	4038917434522
EAN of outer box	4038917434515

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

