



Product description

Order number	70.3040.305
Product description	Pipette tip, 300 µl, Biosphere® plus, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Brand and Finnpipette and products of identical design, 96 pcs./SingleRefill

Product characteristics

Graduation	yes
Volume of work	300 µl
Filter	no
Coding colour	orange
Autoclavable (max. °C)	121 °C

Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	yes
Pyrogen-free	yes
Lot or serial number	Lot number

Packaging

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Pcs. / inner box	480
Pcs. / box	1920
Pcs. / pallet	103680
Depth of case	437 mm
Width of case	351 mm
Height of case	172 mm
Case volume	0,0264 cbm
Weight of case	2,74 kg
EAN of inner box	4038917433136
EAN of outer box	4038917433129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

