



Product description

Order number	70.3050.020
Product description	Pipette tip, 1000 µl, blue, calibration rings, fits Eppendorf, Gilson, Biohit, Socorex, Finnpipette and Brand and products of identical design, 500 pcs./bag

Product characteristics

Graduation	yes
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product blue

Purity & certification

Product in accordance with IVD
MDD or IVDD
Certificate of conformity yes
CE certified CE
Purity standard PCR Performance Tested
Pyrogen-free no
Lot or serial number Lot number

Packaging

Minimum order qty. 5000
Type of inner packaging bag
Pcs. / inner box 500
Pcs. / box 5000
Pcs. / pallet 180000
Depth of case 398 mm
Width of case 398 mm
Height of case 300 mm
Case volume 0,0475 cbm
Weight of case 4,74 kg
EAN of inner box 4038917434584
EAN of outer box 4038917434577

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

