



### **Product description**

Order number	70.3050.110
Product description	Low retention pipette tip, 1000 µl, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Socorex, Finnpipette and Brand and products of identical design, 480 pcs./StackPack

### **Product characteristics**

Graduation	yes
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	1920
Type of inner packaging	StackPack
Pcs. / inner box	1920
Pcs. / box	1920
Pcs. / pallet	149760
Depth of case	303 mm
Width of case	295 mm
Height of case	225 mm
Case volume	0,0201 cbm
Weight of case	2,79 kg
EAN of inner box	4038917424813
EAN of outer box	4038917424806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

