



Product description

Order number	70.3050.205
Product description	Pipette tip, 1000 µl, Biosphere® plus, transparent, calibration rings, fits Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 pcs./box

Product characteristics

Graduation	yes
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	yes
Pyrogen-free	yes
Lot or serial number	Lot number

Packaging

Minimum order qty.	1920
Type of inner packaging	Box
Pcs. / inner box	480
Pcs. / box	1920
Pcs. / pallet	61440
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0,0436 cbm
Weight of case	4,99 kg
EAN of inner box	4038917424752
EAN of outer box	4038917424745

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

