



### **Product description**

Order number	70.3050.375
Product description	Low retention filter tip, 1000 µl, Biosphere® plus, transparent, calibration rings, fits Eppendorf, Gilson and Brand and products of identical design, 96 pcs./SingleRefill

### **Product characteristics**

Graduation	yes
Volume of work	1.000 µl
Filter	yes
Coding colour	blue

### **Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	yes
Pyrogen-free	yes
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Pcs. / inner box	480
Pcs. / box	1920
Pcs. / pallet	69120
Depth of case	512 mm
Width of case	325 mm
Height of case	237 mm
Case volume	0,0394 cbm
Weight of case	4,26 kg
EAN of inner box	4038917424578
EAN of outer box	4038917424561

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

