



Product description

Order number	70.3060
Product description	Pipette tip, 1250 µl, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Socorex, Finnpipette and Brand and products of identical design, 500 pcs./bag

Product characteristics

Graduation	yes
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	no
Lot or serial number	Lot number

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Pcs. / inner box	500
Pcs. / box	5000
Pcs. / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0,0657 cbm
Weight of case	5,72 kg
EAN of inner box	4038917434645
EAN of outer box	4038917434638

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

