



Product description

Order number	70.3060.100
Product description	Pipette tip, 1250 µl, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Socorex, Finnpiquette and Brand and products of identical design, 480 pcs./StackPack

Product characteristics

Graduation	yes
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	no
Lot or serial number	Lot number

Packaging

Minimum order qty.	1920
Type of inner packaging	StackPack
Pcs. / inner box	1920
Pcs. / box	1920
Pcs. / pallet	184320
Depth of case	290 mm
Width of case	267 mm
Height of case	225 mm
Case volume	0,0174 cbm
Weight of case	2,97 kg
EAN of inner box	4038917432207
EAN of outer box	4038917432191

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

