



Product description

Order number	70.3060.110
Product description	Low retention pipette tip, 1250 µl, transparent, calibration rings, fits Eppendorf, Gilson, Biohit, Socorex, Finnpipette and Brand and products of identical design, 480 pcs./StackPack

Product characteristics

Graduation	yes
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product in accordance with MDD or IVDD IVD
Certificate of conformity yes
CE certified CE
Purity standard PCR Performance Tested
Pyrogen-free no
Lot or serial number Lot number

Packaging

Minimum order qty. 1920
Type of inner packaging StackPack
Pcs. / inner box 1920
Pcs. / box 1920
Pcs. / pallet 184320
Depth of case 290 mm
Width of case 267 mm
Height of case 225 mm
Case volume 0,0174 cbm
Weight of case 2,97 kg
EAN of inner box 4038917432238
EAN of outer box 4038917432221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

