



Product description

Order number	70.3060.305
Product description	Pipette tip, 1250 µl, Biosphere® plus, transparent, calibration rings, fits Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 pcs./SingleRefill

Product characteristics

Graduation	yes
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	IVD
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	yes
Pyrogen-free	yes
Lot or serial number	Lot number

Packaging

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Pcs. / inner box	480
Pcs. / box	1920
Pcs. / pallet	69120
Depth of case	512 mm
Width of case	325 mm
Height of case	237 mm
Case volume	0,0394 cbm
Weight of case	4,38 kg
EAN of inner box	4038917432474
EAN of outer box	4038917432467

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

