

MICRO-TOUCH® DENTA-GLOVE® Blue Nitrile

Nitrile Powder-Free Non-Sterile Examination Glove

Issue date 29-02-2020

PRODUCT DESCRIPTION

Material	Nitrile
Colour	Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Fingertip textured
Internal Surface	Polyacrylate coated

PHYSICAL PROPERTIES

Thickness (single)	Finger	0.11
Typical average values (mm)	Palm	0.08
	Cuff	0.07
Minimum length (mm)	240	
Strength (average values)	before ageing	after ageing
Elongation at Break (%)	≥ 700	≥ 650
Force at break (N)	≥ 6N	≥ 6N

PACKAGING AND STORAGE

Packaging	10 x 100 pieces / 1000 pieces per case, 10 x 90 pieces / 900 pieces per case (size:x-large)
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight. Store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCT REFERENCES

Size / product code	x-small 3130050XS	small 31300500S
	medium 31300500M	large 31300500L
	x-large 3130050XL	

KEY FEATURES AND ADVANTAGES

- Better grip due to textured surface
- Comfort and sensitivity
- Easy donning and softer feel
- Minimizes Type I allergy risk
- Specifically designed to meet dentists needs

MANUFACTURING AND SAFETY STANDARDS

AQL (pinholes)	<ul style="list-style-type: none"> • Waterleak tested • Conform European Norm EN 455-1: Inspection Level I AQL 1.5
Protein content	Not applicable: contains no natural rubber latex
Primary skin irritation	Not considered a primary irritant as per FHSA Regulation guidelines 16 CFR 1500
Skin Sensitization	No evidence of delayed dermal contact sensitisation as per ISO 10993-10
Viral Penetration	Passes ISO 16604 using Phi X 174
Sterilization	Not Applicable
CE Mark	Personal Protective Equipment: Category III Medical Device: Class I
Product Standards Conformance	EN 455 parts 1, 2, 3, 4
	EN 420
	EN ISO 374 part 1
	EN 374 part 2
	EN 16523 part 1
	EN 374 part 4
	EN ISO 374 part 5
Manufacturing Standards	
Registering Authority	British Standards Institution (0086): PPE

PICTOGRAMS



The information contained in this document is tentative and might be subject to changes

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu