

# MICRO-TOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup> Blue Nitrile

DENTA-GLOVE<sup>®</sup> Series

## Nitrile, powder-free blue examination glove with ergonomic design for extra comfort

This non-sterile, blue nitrile examination glove delivers the ultimate fit and comfort. Designed with ERGOFORM technology, this glove helps reduce muscle effort that can cause stress on joints, ligaments and tendons. Durable enough for rugged donning while providing great tactile sensitivity and grip for fine, detailed procedures.

### KEY FEATURES AND BENEFITS:

- Innovative film technology provides comfort and fit
- Robust strength helps avoid rips and tears while donning and inner polymer coating facilitates ease of donning
- This blue exam glove offers textured finger surfaces provide a secure grip on instruments in wet or dry conditions
- Tested Safe for Oral Examination, ensuring that the glove is safe for patient contact

### Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Extra comfort during long, sustained procedures
- Protection from Type I latex allergy in HCW's or patients
- Oral examinations
- Standard, moderate risk examination procedures

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup> Blue Nitrile

D E N T A - G L O V E <sup>®</sup> S e r i e s

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.11 / 4.3	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	No
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	DENTA-GLOVE <sup>®</sup> Series

## ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	31300500L	31300500M	31300500S	3130050XL	3130050XS

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *90 gloves per dispenser (Size XL) *900 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

Category III, EN 420:2003 + A1:2009, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 9001, ISO 13485

## Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com