



### Product description

Order number	62.612
Product description	Tube 5 ml, 92 x 15.3 mm, PP, round base (open), with natural false bottom and assembled screw cap, with print, 100 pieces packed in the bag

### Product characteristics

Label/ Print	With print
Colour of print/label	white
Graduation	yes
Base shape	Flat false bottom

### Size

Sample volume	5 ml
Diameter	15,3 mm
Height of product	92 mm
Length of product	97 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of case	5,69 kg
EAN of inner box	4038917082884
EAN of outer box	4038917001908

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

