



## MICRO-TOUCH® DENTA-GLOVE® HYDRACARE®

### SKIN CARE COMBINES WITH PROTECTION FOR DENTAL SPECIALISTS

Developing best-in-class hand protection begins with our underlying focus on users – what they do, how they do it and how we can help them do it better. This approach enable us to purpose-build our gloves, going beyond safety and protection to promote complete hand health.

Problems facing the dental community: Hand fatigue is a serious issue.

Prolonged muscle effort over long period of times can cause hand fatigue, which has been associated with Carpal Tunnel Syndrome, as well as other hand injuries. These conditions are serious, painful and can reduce productivity, especially for dentists and hygienists whose job makes wearing disposable exam gloves an essential requirement.

To minimise the effects of continuous glove use on hands, Ansell offers MICRO-TOUCH® Denta-Glove® HydraCare®: this powder-free latex glove contains an interior glycerol moisturiser that helps maintain skin suppleness. This has a proven corrective effect on defective skin, restoring hydration and elasticity.

Single chlorination combined with silicone and the moisturiser guarantee easy donning, even with damp hands. A pale green colour ensures easy product identification.

### PRODUCT DESCRIPTION

Material	Natural Rubber Latex
Colour	Pale Green
Shape	Ambidextrous
Cuff	Beaded
External Surface	Micro-Textured
Internal Surface	Slightly chlorinated, Treated with silicone and moisturizer

### SIZE / PRODUCT CODE

<b>x-small</b>	553621	<b>small</b>	553622
<b>medium</b>	553623	<b>large</b>	553624
<b>x-large</b>	553625		

MICRO-TOUCH® gloves are designed and engineered for a broad range of hospital staff, emergency medical technicians, healthcare facility personnel, and private practitioners seeking the protection provided by examination gloves for the variety of settings where they face challenging risks every day. MICRO-TOUCH® gloves offer hand protection to give users the confidence and professional freedom to do their jobs.

Natural Latex Powder-Free Non-Sterile Examination Glove

MD Class I

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at [www.ansell.eu](http://www.ansell.eu).

**Ansell Healthcare Europe N.V. (European Head Office)**

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail [info@ansell.eu](mailto:info@ansell.eu)