



### **Product description**

Order number	60.611.010
Product description	Tube 5 ml, 92 x 15.3 mm, with false bottom, flat, PP, without cap

### **Product characteristics**

Graduation	no
Base shape	Flat false bottom

### **Size**

Sample volume	5 ml
Diameter	15,3 mm
Height of product	92 mm
Length of product	92 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of case	4,68 kg
EAN of inner box	4038917114820
EAN of outer box	4038917014427

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

