

Medical face masks - Type IIR - High splash protection

M E D I C A L F A C E M A S K S

Medical face masks help prevent large particles expelled by the wearer (e.g. spit, mucous) from reaching the patient or work environment.

Protecting healthcare staff and patients from infectious diseases

Medical face mask, Type IIR - EN14683, features include:

- Pleat style with ear loops or ties
- Protective four-layer protection
- Splash resistant layer protection against blood and bodily fluids
- Available in a variety of colours and styles

Medical face masks help prevent large particles expelled by the wearer (e.g. spit, mucous) from reaching the patient or work environment and help protect the wearer from exposure to blood and/or bodily fluids from the patient or work environment.

Type IIR - EN14683 masks are appropriate for situations where exposure to blood and/or bodily fluids from the patient or work environment is at risk for the healthcare worker.

KEY FEATURES & BENEFITS

- Pleat style with ear loops or ties
- Protective four-layer protection
- Splash resistant layer protection against blood and bodily fluids
- Available in a variety of colours and styles

RECOMMENDED FOR

- Enhancing infection control
- Preventing the risk of cross-contamination



SANDEL[®] Anticipate the Unexpected

Innovative SANDEL[®] Safety Solutions are developed through collaboration with clinical professionals to help prevent injuries, infections, and errors in high-risk healthcare environments. SANDEL products provide the confidence that healthcare professionals need to perform at their best and protect patients from preventable harm.

Medical face masks - Type IIR - High splash protection

M E D I C A L F A C E M A S K S

PRODUCT DESCRIPTION

Color	Light Blue, Green
Shape	Flat Pleated
String Type	Ear loops, Ties

MANUFACTURING AND SAFETY STANDARDS

Lot Traceability	Yes
Reusable	No, Disposable Product, Single Use Only
Resterilization	Not Applicable
Shape	Flat Pleated
String Type	Ear loops, Ties
Standards	EN 14683:2014, EN ISO 13485, ISO 9001, ISO 22609:2004

PACKAGING AND STORAGE

Packaging Overview	
Case Quantity	400, 400, 400, 400, 100, 100
Shelf Life	5 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition

Contact your Ansell representative for ordering or more information.

North America

Tel: (866) 764-3327
 Email: sandel-cs@ansell.com
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Latin America & Caribbean

Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
 Brazil Tel: +55-11-3356-3100
 Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa

Tel: +32 (0) 2 528 74 00
 Email: info@anselleu

Asia Pacific

Tel: +603 8310 6688
 Email: apac.medical@ansell.com

Australia & New Zealand

Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com