



Product Catalogue



UltrageL was established in 1997, has created this company with the purpose of producing – in Hungary for the following Medical Accessories

- ◆ Contact Gels for Ultrasound, IPL, ECG
- ◆ Sports Gels
- ◆ Hand Sanitizer Gels
- ◆ Disposable and Reusable ECG electrodes
- ◆ Defibrillator Electrodes
- ◆ ECG Patient Cables
- ◆ Medical Chart Papers
- ◆ Videoprinter Papers
- ◆ Disposable Viral and Bacterial Filter for Pulmonary Function Test

The ultrasound- and ecg-gels produced and marketed by us on the basis of our own formula and the individual production technology we developed are well known and widely distributed all over in the international market.

The excellent quality of our gels, and other medical accessories is acknowledged in many countries around the world, since UltrageL's products are present in Asia, Africa, Europe, Far & Middle Eastern Countries as well. Thanks to the commercial structure of the company covers the entire national territory in a capillary manner and it exports to more than 34 countries worldwide.

Proper product quality of the produced medical items is ensured through applying precise production technology and strict controls, based on a closed production system and on a continuously developing equipment park. Controlling under the International System Standards the digital PH-measuring, computerised viscosity control, digital water quality examination and microbiological report also assure best production level.

In compliance with our traditions our goal is to provide our customers top quality services, also meeting the high quality standards. UltrageL is a Company with Certified Quality Management System in accordance with ISO 9001:2008 quality assurance system from TÜV Rheinland, EN ISO 13485:2012. All UltrageL products are CE marked according to Directive 93/42/EEC.

The production takes places on 2000 square meters in Budapest, where is together international and domestic sales. We are also taking care of R&D, warehousing, logistics, and administration purposes.



HiGel Hand Sanitizer

70% alcoholic content
Laboratory tests prove its effectiveness.

Bactericid, MRSA, fungicid, virucid, tuberculocid.


Practical:

- ◆ in households,
- ◆ hospitals,
- ◆ retirement homes,
- ◆ nursery schools,
- ◆ sports complexes,
- ◆ private practices
- ◆ and in the everyday when hand sanitizing urgently essential.

Leave hands feeling refreshed with pleasant lemon fragrance and without sticky residue. Thanks to its added components keeps the skin delicate.

Minimal effect time: 30 sec
Virus inactivation time: approx. 2 minutes

Versions:100 ml, 300 ml, or according to request.



AquaUltra Basic


Traditional quality gel with blue colour, high viscosity (82000 RPM5-STD 18o), easily dissolving contact material. Recommended for all procedures.



AquaUltra Steril Gel

Sterilized ultrasound transmission gel is recommended when sterility is indispensable (e.g.: diagnostic ultrasound examination, in operating room, catheterization)

Available in 5g, 20 g pouches




Aqua Laser

Specifically high viscosity gel, water-soluable, colourless and odorless contact material, ideal for Cosmetic Laser, IPL, Therapy.



Transept Spray

Manufacturers recommend to use it after each examination. It cleans the probe without scratching and scraping it's surface.




AquaUltra Clear

Colourless and fragrance free. Easily dissolving contact material with high viscosity level.



AquaUltra Gels for Ultrasound

- Perfect conductivity
- Hypoallergenic
- Non-irritating
- Free from salt and formaldehyde
- Non-mucous, non-adhesive,
- Water-soluble
- Only food colour pigments used
- Rapid transmission speed



Conti ECG Gel

Specifically high-conductivity contact gel for all electromedical procedures. (ECG, EMG, defibrillation, and biofeedback).



Accessories for rigid tank

Eco Dispenser: lower capacity pump for canister.

Max Dispenser: high-capacity pump recommended for continuous use.



Conti Cream

Specifically high-conductivity contact cream for all electromedical procedures. (Special for EEG, EMG).



Conti Spray

Specifically high electric conductivity electrolyte spray. Hypoallergenic, antibacterial, economical and easy to spread. After using, it can easily be wiped off. Ideal for stress test systems.



Conti Gels for Electrocardiographic Diagnostics

- High electric conductivity and quality.
- Low electric resistance
- Recommended for short- and long term recordings
- Bacteriostatic and non-irritating (pH from: 6,5 to 7)
- Water-soluble, not greasy
- Hypoallergenic
- Optimal viscosity
- Easy to wipe off




Under the most difficult conditions a maximum result is guaranteed in stability and precision as well. The electrodes line is totally complete and therefore meets any monitoring requirements. Both for adult/pediatric use, for normal and stress tests. Base made of **Foam** and **Micropore**, in the form of oval, rectangular and round shaped. With all these sortiment we provide excellent complementary or similar articles. As a result of the continuous research and information exchange, with the end-customer, we are able to manufacture high quality and reliable products.




Adapters for disposable electrodes

Snap adapters between clip electrodes and cables with 4 mm plug in 5 colours. Alligator adapters with socket connector. Universal adapters with 4 mm socket for snap & alligator connection in 5 colours.



Disposable ECG electrodes

Foam based, Ag/AgCl connector for normal monitoring, pediatric use, stress tests and holter. Diameter from 25 mm to 60 mm. Radio-lucent electrodes with carbon fibres sensor. Round, oval, rectangular shapes.



	Solid gel	Liquid gel	Monitoring	Stress test	Dual
F 25	•		•		
F 50	•		•		
F 55	•		•		
FS 50	•	•	•	•	
F 5035	•		•	•	•
F 37x41	•		•		
FS 37x41	•			•	
F 32x40		•	•		
F 32x40	•		•		
FS 32x40	•			•	

Electrodes & Accessories

Suction chest & Clamp electrodes, Rubber belts

Extremely high conductivity thanks to the Ag/AgCl sensors which guarantees the highest reliability in signal transmission. Universal connection. Available in two different sizes: 15 mm, 24 mm. Packaging: bag 6 pcs, 4 pcs.



Patient cables

Large range of 4 mm cables with plug or snap connection. The available connectors are compatible with the most ECG devices, presented on the market. For each model, with or without protection from defibrillation discharges cables are available. 3, 5 or 10 leads, service- and detachable cables.

Compatible with the following devices:

Esaote, ET & Cardioline, Aspel, Cardiorapid, Hellige, Siemens/Burdick, Kontron, Schiller, Fukuda Denshi, BTL, Innomed, Medcor, Carditek, Dr Lee, Progetti, Macquare, Hewlett Packard, Bosch, Mortara, Picker, Cambridge, Kenz, Marquette, etc.



Defi pads


Disposable multifunction electrodes for defibrillation, synchronized cardioversion, transcutaneous cardiac stimulation and ECG monitoring, with direct connection to cables and defibrillators. Every model is equipped with a cable of the right length and safety antishock connector. Thanks to its technical features they grant a low impedance and stimulation threshold, with a very good patient compliance.

Primary packaging: 1 pair in polyethylene and aluminium sealed pouch

Sale packaging: 10 pairs in carton box
Pediatric version is also available.

Compatible with the following devices:

Philips, Zoll, Medtronic – Osato Bexen – Cardioline, Agilent – Laerdal Medical, Schiller – Esaote, Nihon Kohden – Welch Allyn – Innomed, Shiller, CU Medical Systems, Metrax, Primedic, GE – Esaote, Cmos Drake, Nihon Kohden



Recording papers for diagnostic

Manufacturing all types of medical recording charts in roll and Z-fold using the latest technology printing machines, resulting products with absolute precision. From the single channel machines to the most complex appliances for stress test or particular monitoring of coronary units.

We guarantee the absolute reliability of our recording papers with respect to all the instruments available commercially, even the less common. At Ultrigel we have a wide range of anti fade thermal papers with guaranteed image life span of up to 20 years to suit every need in the medical sector. Ultrigel is dedicated to manufacturing any required medical recording charts to match all types of equipment, used in hospitals, laboratories and a wide range for industries.


Wide range of the papers for various of sectors, including:

- ♦ electrocardiographs
- ♦ intensive care equipments
- ♦ electroencephalographs
- ♦ neurology recorders
- ♦ foetal monitors
- ♦ gynecology recorders
- ♦ spirometrical and audiometrical equipments
- ♦ laboratory analysers



Wide range of videoprinter papers and continuous availability of original materials. Original & compatible Sony and Mitsubishi videoprinter papers

- ♦ black & white
- ♦ colours



Recording papers for diagnostic & videoprinter papers

Electrodes & Accessories



Deep Ice Cooling Sport Gel

This new gel provides controlled cooling without freezing or irritation and has exceptional glide for point therapy.

Effective for acute conditions or injuries that have just occurred, in case of inflammation of joints and muscles.

Relieves minor aches, simple backaches, strains, bruises and sprains. In case of muscle soreness, it's relaxant.

Natural components are menthol and camphor.

No artificial colours. Produced without parabens and SLS.

Available sizes: 260 g, 1000 g

HOT & COLD Gelpack

Hot/cold painkiller treatment. It is available in two different sizes: 13,5x18 cm, 13,5x28 cm, or according to request.

Deep Training Thermo Sport Gel

Ultragel's unique formula development blending 4 natural plant extracts – rosemary, eucalyptus, lemon and cinnamon – with safe active ingredients.

The result is an elegant gel which is also able to warm up the muscles before training. Reduces the warm up time and the risk of injuries.

Effects are immediate and will not stain clothing.

No artificial colours. Produced without parabens and SLS.

Available sizes: 260 g, 1000 g, or according to request.

Sport Gels



ULTRAGEL HUNGARY 2000 Kft.

Hungary, 1211 Budapest, Tekercselő u. 12. P.: +36 1 278 3057 F: +36 1 278 3055 E: info@ultragel.hu
www.ultragel.hu