

PROFEEL[®]

DHD[™] Platinum

Powder Free Surgical Gloves



Damp Hand Donning - DHD[™] • Natural Rubber Latex • Powder Free

www.wrpworld.com

CE
2797
BSI (2797)
1066 EP, NL

EN ISO 374-1:2016
Type B

KPT

EN ISO 374-5:2016

VIRUS

PROFEEL®

DHD™ Platinum

Powder Free Surgical Gloves

- Natural Rubber Latex's intrinsic softness provides an exceptionally comfortable and flexible glove.
- The glove has been specially formulated to enable effortless double gloving.
- Damp Hand Donning (DHD)- The inner surface of the glove is coated with a proprietary polymer to enhance damp hand donning
- A non-slick skin-grip surface finish provides excellent instrument grip, under both wet and dry conditions.
- A superior tactile sensitivity provides optimum 'feel'.
- Specially designed high tack band near the cuff of the glove to securely hold onto the wearer's gown to minimize roll down.
- A reinforced beaded cuff feature ensures easy donning while preventing tearing.

PRODUCT INFORMATION

Article No.	: WRP SG NR PF 125
Style No.	: WRP SG NR PF 125 - PTN DHD
Material	: Natural Rubber Latex
Coating	: Damp Hand Donning (DHD)
Weight per pair, size 7.5, g	: 23.4 ± 0.6
Sizes	: 5.5 - 9.0
Palm Thickness, mm (min)	: 0.17
Length (mm)	: 290 ± 10
Colours	: WHITE

Glove Size	Palm Width (mm)
5 ½	72 ± 4
6	77 ± 5
6 ½	83 ± 5
7	89 ± 5
7 ½	95 ± 5
8	102 ± 6
8 ½	108 ± 6
9	114 ± 6

PRODUCT CONFORMANCE:

Medical Device: In compliance with European Medical Device Directive 93/42/EEC (Class IIa), EN 455: Part 1/2/3/4, EU Regulation 2016/425, EN 420:2003 +A1 2009, EN ISO 374-1:2016, EN ISO 374-2:2014, EN ISO 374-4:2013, EN 374-5:2016, EN 16523-1:2015, ASTM D3577, FDA 510 K

QUALITY ASSURANCE:

Manufacturing process is in compliance with US FDA Quality System Regulation (QSR), ISO 9001 Quality Management System and BS EN ISO 13485 Quality Systems



 Your Partner In Protection™

EC REP

REMESCO Handelsges.m.b.H.

Grinzinger Allee 5/25, 1190 Vienna, AUSTRIA

Office : +43-1-328-50-88

Fax : +43-1-328-29-49

Email : office@remesco.com

WRP Asia Pacific Sdn Bhd 147817V

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi, 43900 Sepang, Selangor Darul Ehsan, MALAYSIA

Telephone : +60-3-8706 1486

Fax : +60-3-8706 1557

Email : customer_wrp@wrpworld.com
sales@wrpworld.com

Website : www.wrpworld.com